



Digital og automatisk urinmåling for alle sykehusavdelinger



Sippi™ – Digital og automatisk urinmåling

Brukervennlig og tidsbesparende urinmålere til alle sykehusets avdelinger.

Produktinformasjon

Baseenhet

- Digital og automatisk urinmåling
- Måler både døgndiurese og timediurese
- Berøringsfri kapasitiv måleteknologi
- Mulighet for trådløs dataoverføring til sykehusets PDMS-system
- Måleintervall 0–2000 ml/t
- Forventet livssyklus er 3 år
- 3 stk. AA-batterier (batteritid 4 måneder)
- Gjeldende væskedøgns diurese er alltid tilgjengelig for kontroll og oppfølging

Urinpose

- Totalvolum 2000 ml
- Kan tømmes
- Nålefri prøvetakingsport
- Anbefalt brukstid opptil 7 døgner
- Steril (EO)
- Kapsel med silikonolje er plassert i målekammeret (SippCoat™). Ved kontakt med urin frigjøres silikonoljen, som motvirker biofilm.^{1,2}



Ordreinformasjon

Art.nr.	Beskrivelse	Pakning
1104	Sippi baseenhet BLE – digital urinmåler	1 stk
1203	Sippi urinposer	20 stk

Overholder regulerende krav i:
LVFS 2003:11, MDD93/42/EEG

Referanser: 1. Slettengren M, et al. "Making medical devices safer: impact of plastic and silicone oil on microbial biofilm formation." *The Journal of Hospital Infection*. Vol. 106,1 (2020): 155-162. 2. Slettengren M, et al. "Silicone oil decreases biofilm formation in a capacitance-based automatic urine measurement system." *Sensors*. 2021; in press