



Digital och automatisk urinmätning för sjukhusets alla avdelningar



Sippi™ – Digital och automatisk urinmätning

Användarvänlig och tidsbesparande urinmätare för sjukhusets alla avdelningar.

Produktinformation

Basenhet

- Digital och automatisk urinmätning
- Mäter både dygnsdiures och timdiures
- Beröringsfri kapacitiv mätteknologi
- Möjlighet för trådlös dataöverföring till sjukhusets PDMS
- Mätintervall 0–2000 ml/h
- Förväntad livscykel är 3 år
- 3 st AA-batterier (batteritid 4 månader)
- Nuvarande vätskedygns diures alltid tillgänglig för kontroll och uppföljning

Urinpåse

- Totalvolym 2000 ml
- Tömningsbar
- Nålfri provtagningsport
- Rekommenderad användningstid upp till 7 dygn
- Steril (EO)
- Kapsel innehållande silikonolja är placerad i mätkammaren (SippCoat™). Vid kontakt med urin frigörs silikonoljan, vilket motverkar biofilm.^{1,2}



Orderinformation

Art.nr.	Beskrivning	Förp
1104	Sippi Base Unit BLE – digital urinmätare	1 st
1203	Sippi urinpåsar	20 st

Möter de regulatoriska kraven i:
LVFS 2003:11, MDD93/42/EEG

References: 1. Slettengren M, et al. "Making medical devices safer: impact of plastic and silicone oil on microbial biofilm formation." *The Journal of Hospital Infection*. Vol. 106,1 (2020): 155-162. 2. Slettengren M, et al. "Silicone oil decreases biofilm formation in a capacitance-based automatic urine measurement system." *Sensors*. 2021; in press